

M-TURBO®

DESIGNED FOR THE POINT OF CARE



THE M-TURBO – FOR HIGH RESOLUTION IMAGING ACROSS THE POINT OF CARE

Engineered for striking image quality, durability and ease of use, our flagship product, the M-Turbo® ultrasound system, has earned customer satisfaction ratings of 99%*. Proprietary algorithms optimise multiple imaging parameters so you get essential information quickly and easily. The M-Turbo ultrasound system offers an advanced set of features with a wide array of connectivity options that seamlessly connects you to hospital information networks and your own PC.

*Verispan surveys 2008

TURBO TECHNOLOGY FROM SONOSITE

The M-Turbo system's 16-fold increase in processing power generates dramatic improvements in image quality by simultaneously running multiple advanced algorithms:

SonoADAPT™ Tissue Optimisation – eliminating complicated manipulation of multiple controls.

SonoHD™ Imaging Technology – reducing speckle noise and other image artifacts while preserving and sharpening tissue information.

SonoMB® Multi-beam Imaging – increasing resolution of small structures and enhancing border delineation.

Advanced Needle Visualisation – aiding needle visualisation while maintaining striking image quality of the target and surrounding anatomy.

ColorHD™ Technology – increasing colour performance, sensitivity and frame rates for more diagnostic information.


POINT-OF-CARE APPLICATIONS INCLUDE:

Anaesthesia, Critical Care, Cardiology, Cardiovascular Disease Management, Emergency Medicine, Musculoskeletal, Ob/Gyn, Radiology, Vascular, Surgery, Shared Service, Women's Health

To learn more about our products, contact us at **01462 444800** or visit www.sonosite.com/products/m-turbo



Advanced Needle Visualisation – Saphenous Nerve Block, using the HFL50x/15-6 MHz transducer



Carotid Artery, using the HFL38x/13-6 MHz transducer



M-TURBO®

DESIGNED FOR THE POINT OF CARE



Capture quality video clips up to 60 seconds long

Battery powered and wireless capability for true mobility

PC- and Mac-friendly for effortless data management with 2 high-speed USB 2.0 ports

Weighs 3.4 kg/7.5 lbs. with battery

Backlit keyboard is easier on the eyes

Magnesium case for lightweight strength

Increased application flexibility with optional Triple Transducer Connect



MAIN FEATURES AT A GLANCE

- Premium image quality
- Drop tested at 91.4 cm/3 feet
- Splash resistant user interface
- Quick boot-up time
- Easy to operate
- Wireless connectivity
- No annual service contract



SonoSite Ltd
European Headquarters
Alexander House, 40A Wilbury Way, Hitchin, Herts SG4 0AP, United Kingdom
Tel: +44 1462-444800 Fax: +44 1462-444801 E-mail: europe@sonosite.com

www.sonosite.com/products/m-turbo

SonoSite Worldwide Offices

SonoSite Australasia Pty Ltd – Australia – New Zealand	1300-663-516 0800-888-204 +55 61 8118 7100 +1 888-554-5502
SonoSite – Brazil	+49 91-31-97-79-29-0
SonoSite Canada Inc. – Canada	+34 91-640-49-11
SonoSite GmbH – Germany	+1 425 951-1200
SonoSite Iberica SL – Spain	+91 124-288-1100
SonoSite, Inc. – USA	+39 02-98490885
SonoSite India Pvt Ltd – India	+81 3-5304-5337
SonoSite Italy S.r.l. – Italy	+86 21-5239-6693
SonoSite KK – Japan	+33 1-69-18-96-30
SonoSite Ltd – China – Shanghai	
SonoSite SARL – France	

M-Turbo, SonoSite, the SonoSite logo, and other trademarks not owned by third parties are registered or unregistered intellectual property of SonoSite, Inc., a subsidiary of FUJIFILM.

Mac is a trademark of Apple Inc., registered in the U.S. and other countries.

M-TURBO®

SYSTEM SPECIFICATIONS



THE M-TURBO® ULTRASOUND SYSTEM FROM SONOSITE



SYSTEM SPECIFICATIONS:

System weight	3.04 kg/6.7 lbs.
Dimensions	30.2 cm x 27.4 cm x 7.9 cm /11.9" x 10.8" x 3.1" (L x W x H)
Display	26.4 cm/10.4" diagonal LCD (NTSC or PAL)
Architecture	All-digital broadband
Dynamic range	Up to 165 dB
Gray scale	256 shades
HIPAA compliance	Comprehensive tool set

IMAGING MODES AND PROCESSING:

Broadband, Multifrequency Imaging:

- 2D / Tissue Harmonic Imaging / M-Mode
- Velocity Colour Doppler / Colour Power Doppler
- PW, PW Tissue Doppler and CW
- Doppler angle, correct after freeze

Image processing:

- SonoADAPT™ Tissue Optimisation
- SonoHD™ Imaging Technology
- Advanced Needle Visualisation (SonoMBe™ Imaging)
- Dual Imaging, Duplex Imaging, 2x pan/zoom capability, Dynamic range and gain

USER INTERFACE AND REMAPPABLE CONTROLS:

- Softkeys to drive advanced features
- Programmable A and B keys: each can be assigned by the user for increased ease of use
- Alphanumeric elastomeric QWERTY keyboard
- Track pad with select key for easy operation and navigation
- Doppler controls: angle, steer, scale, baseline, gain and volume
- Image acquisition keys: review, report, Clip Store, DVD, save
- Dedicated AutoGain and exam keys to allow quick activation

TRANSDUCERS:

Broadband and Multifrequency:

- Linear Array, Curved Array, Phased Array, Multiplane TEE and Micro-Convex

Single Frequency:

- Cardiac Static Pencil

APPLICATION SPECIFIC CALCULATIONS:

OB/GYN/Fertility: Diameter/ellipse measurements, volume, ten follicle measurements, estimated foetal weight, established due date, gestational age, last menstrual period, growth charts, user-defined tables, multiple user-selectable authors, ratios, amniotic fluid index, patient report, humerus and tibia measurement and charts

Vascular: Diameter/ellipse/trace measurements, volume, volume flow, percent diameter and area reduction, Lt/Rt CCA, ICA, ECA, ICA/CCA ratio, time average mean (TAM), peak trace, ICA/CCA ratio, angle correction, patient report

CIMT (Carotid Intima Media Thickness):

Embedded SonoCalc® IMT software (optional)
- automatic edge detection with mean and maximum thickness reporting

Cardiac: Automated Cardiac Output package and patient report including: ventricular, aortic and atrial measurements; ejection fraction, volume measurements, Simpson's rule, continuity equation, pressure half-time and cardiac output; PA AT, TV E, A, PHT, TVI, MV time, Pulm Veins

Transcranial Doppler (TCD): Complete TCD package including Time Average Peak (TAP)



M-TURBO®

SYSTEM SPECIFICATIONS



ONBOARD IMAGE AND CLIP STORAGE/REVIEW:

- 8 GB internal Flash memory storage capability
Potential to store 30,000 images or 960 2-second clips
- Clip Store capability (maximum single clip length: 60 seconds)
- Clip Store capability via either number of heart cycles (using the ECG) or time base
Maximum storage in ECG beats mode is 10 heart cycles. Maximum storage in time base mode is 60 seconds
- Cine review up to 255 frame-by-frame images

POWER SUPPLY:

- System operates via battery or AC power
- Rechargeable lithium-ion battery
- AC: universal power adapter, 100-240 VAC, 50/60 Hz input, 15 VDC output

MEASUREMENT TOOLS, PICTOGRAMS AND ANNOTATIONS:

2D: Distance calipers, ellipse and manual trace
Doppler: Velocity measurements, pressure half time, auto and manual trace

M-Mode: Distance and time measurements, heart rate calculation

User-selectable text and pictograms

User-defined, application-specific annotations
Biopsy guidelines

EXTERNAL DATA MANAGEMENT AND WIRELESS:

DICOM® Image Management (TCP/IP):

- Print and Store, Modality Work List

Storage Commit:

- Modality, Perform, Procedure Step

PC Workstation Image Management (TCP/IP, USB):

- SiteLink™ Image Manager – allows transfer, archiving, viewing and printing of high resolution bitmap images/clips, and batch compression to JPEG on PCs
- Direct writing capability to USB 2.0 mass storage removable media (PC and MAC compatible)
- Supported export formats are: MPEG-4 (H.264), JPEG, BMP, and HTML
- SonoSite® Education Key™ training video compatible

SonoSite Workflow Solutions (SWS™):

- A suite of specialised tools to streamline exam management for billing, credentialing, EMR and image archiving

EXTERNAL VIDEO AND AUDIO:

- S-video (in/out) to VCR or DVD for record and playback
- RGB or DVI output to external LCD display
- Composite video output (NTSC/PAL) to VCR or DVD, video printer or external LCD display
- Audio output
- Integrated speakers

H-UNIVERSAL™ STAND AND PERIPHERALS:

- Transducer and gel holders
- Optional Triple Transducer Connect (TTC) to quickly activate transducers electronically
- Optional foot switch

OPTIONAL PERIPHERALS:

Printers: Medical-grade black and white or colour

External storage devices: Medical-grade DVD

External data input devices: Bar code reader

ECG module: 3-lead ECG – works with standard ECG leads and electrodes

External analog ECG input also available

USB bar code reader

M-Turbo, SonoSite, the SonoSite logo, and other trademarks not owned by third parties are registered or unregistered intellectual property of SonoSite, Inc., a subsidiary of FUJIFILM.

Mac is a trademark of Apple Inc., registered in the U.S. and other countries.

DICOM is the registered trademark of the National Electrical Manufacturers Association for its standards publications relating to digital communications of medical information.



SonoSite Ltd

European Headquarters
Alexander House, 40A Wilbury Way,
Hitchin, Herts SG4 0AP, United Kingdom
Tel: +44 1462-444800 Fax: +44 1462-444801
E-mail: europe@sonosite.com

www.sonosite.com/products/m-turbo

SonoSite Worldwide Offices

SonoSite Australasia Pty Ltd – Australia	1300-663-516
– New Zealand	0800-888-204
SonoSite – Brazil	+55 61 8118 7100
SonoSite Canada Inc. – Canada	+1 888-554-5502
SonoSite GmbH – Germany	+49 91-31-97-79-29-0
SonoSite Iberica SL – Spain	+34 91-640-49-11
SonoSite, Inc. – USA	+1 425 951-1200
SonoSite India Pvt Ltd – India	+91 124-288-1100
SonoSite Italy S.r.l. – Italy	+39 02-98490885
SonoSite KK – Japan	+81 3-5304-5337
SonoSite Ltd – China – Shanghai	+86 21-5239-6693
SonoSite SARL – France	+33 1-69-18-96-30

M-TURBO®

TRANSDUCERS
















READY FOR A COMPLETE RANGE OF APPLICATIONS



SonoSite designs, manufactures and tests transducers in-house with real-world customer needs in mind. Our transducers exceed stringent military specifications for drop testing and splash resistance, so you can use them with confidence in the most demanding of environments. Maybe that's why SonoSite is the only ultrasound company to offer an industry leading five-year warranty on the transducers it manufactures.



<p>L38xi</p>  <p>10-5MHz Linear</p> <p>Applications:</p> <ul style="list-style-type: none"> breast CIMT musculoskeletal nerve small parts vascular venous <p>Scan depth:</p> <p>9 cm</p>	<p>HFL38x</p>  <p>13-6MHz Linear</p> <p>Applications:</p> <ul style="list-style-type: none"> breast CIMT musculoskeletal nerve small parts vascular venous <p>Scan depth:</p> <p>6 cm</p>	<p>HFL50x</p>  <p>15-6MHz Linear</p> <p>Applications:</p> <ul style="list-style-type: none"> breast musculoskeletal nerve small parts <p>Scan depth:</p> <p>6 cm</p>	<p>L25x</p>  <p>13-6MHz Linear</p> <p>Applications:</p> <ul style="list-style-type: none"> musculoskeletal nerve superficial vascular venous ophthalmic <p>Scan depth:</p> <p>6 cm</p>	<p>C11x</p>  <p>8-5MHz Curved</p> <p>Applications:</p> <ul style="list-style-type: none"> abdominal neonatal nerve vascular cardiology (Vet) <p>Scan depth:</p> <p>10 cm</p>	<p>C60x</p>  <p>5-2MHz Curved</p> <p>Applications:</p> <ul style="list-style-type: none"> abdominal musculoskeletal nerve ob gyn <p>Scan depth:</p> <p>30 cm</p>
<p>ICTx</p>  <p>8-5MHz Curved</p> <p>Applications:</p> <ul style="list-style-type: none"> ob gyn <p>Scan depth:</p> <p>13 cm</p>	<p>P21x</p>  <p>5-1MHz Phased</p> <p>Applications:</p> <ul style="list-style-type: none"> abdominal cardiology ob orbital TCD <p>Scan depth:</p> <p>35 cm</p>	<p>P10x</p>  <p>8-4MHz Phased</p> <p>Applications:</p> <ul style="list-style-type: none"> paed. abdominal paed. cardiology neonatal head <p>Scan depth:</p> <p>14 cm</p>	<p>SLAx</p>  <p>13-6MHz Linear</p> <p>Applications:</p> <ul style="list-style-type: none"> musculoskeletal nerve superficial vascular venous <p>Scan depth:</p> <p>6 cm</p>	<p>TOEx/TEEx</p>  <p>8-3MHz Multi</p> <p>Applications:</p> <ul style="list-style-type: none"> adult cardiology multiplane transoesophageal 180° rotation of the imaging plane, providing a 360° field of view <p>Scan depth:</p> <p>18 cm</p>	<p>D2x</p>  <p>2MHz Pencil</p> <p>Applications:</p> <ul style="list-style-type: none"> cardiology (CW) <p>Scan depth:</p> <p>N/A</p>
<p>L52x (Vet)</p>  <p>10-5MHz Linear</p> <p>Applications:</p> <ul style="list-style-type: none"> musculoskeletal ob vascular <p>Scan depth:</p> <p>15 cm</p>	<p>Needle guide kits available with the following transducers – L38xi, HFL38x, HFL50x, C60x, ICTx, P10x and P21x. A transverse needle guide is available with the L25x transducer.</p> <p>For detailed information, please visit www.sonosite.com/products/m-turbo or speak with your SonoSite customer representative.</p> <p>M-Turbo, SonoSite, the SonoSite logo, and other trademarks not owned by third parties are registered or unregistered intellectual property of SonoSite, Inc., a subsidiary of FUJIFILM.</p> <div style="display: flex; justify-content: flex-end; align-items: center;"> <div style="background-color: #76b82a; color: white; padding: 5px; margin-right: 10px;"> 5 YEAR </div> <div style="background-color: #34495e; color: white; padding: 5px;"> TECHNOLOGY DRIVEN 5-YEAR STANDARD WARRANTY MADE IN THE USA </div> </div>				



SonoSite Ltd
European Headquarters
Alexander House, 40A Wilbury Way, Hitchin, Herts SG4 0AP, United Kingdom
Tel: +44 1462-444800 Fax: +44 1462-444801 E-mail: europe@sonosite.com

www.sonosite.com/products/m-turbo

SonoSite Worldwide Offices
SonoSite Australasia Pty Ltd – Australia
– New Zealand

SonoSite – Brazil
SonoSite Canada Inc. – Canada
SonoSite GmbH – Germany
SonoSite Iberica SL – Spain
SonoSite, Inc. – USA
SonoSite India Pvt Ltd – India
SonoSite Italy S.r.l. – Italy
SonoSite KK – Japan
SonoSite Ltd – China – Shanghai
SonoSite SARL – France

1300-663-516
0800-888-204
+55 61 8118 7100
+1 888-554-5502
+49 91-31-97-79-29-0
+34 91-640-49-11
+1 425 951-1200
+91 124-288-1100
+39 02-98490885
+81 3-5304-5337
+86 21-5239-6693
+33 1-69-18-96-30