

# MX8 Cardiovascular Ultrasound System

## Performance Specifications

### Application

Pediatric  
Cardiac  
Vascular

### Transducer Types

Linear array  
Phased array

### Imaging Modes

B-Mode  
THI and PSH™ (Phase Shift Harmonic Imaging)  
M-Mode/Color M-mode  
Free Xros M™ (Anatomical M-mode)  
Free Xros CM™ (Curved Anatomical M-mode)  
Color Doppler Imaging  
Power Doppler Imaging/Directional PDI  
Pulsed Wave Doppler  
Continuous Wave Doppler  
TDI  
Low MI Contrast (Myocardium Contrast Imaging)  
Tissue Tracking QA  
Stress Echo  
Strain Elastography

### Standard Features

B-Mode  
THI and PSH™  
M-Mode  
Color Doppler Imaging  
Power Doppler Imaging and Directional PDI  
Pulsed Wave Doppler  
iBeam™ (Spatial Compound Imaging)  
iClear™ (Speckle Suppression Imaging)  
iTouch™ (Auto Image Optimization)  
Echo Boost™  
Zoom/iZoom (Full Screen Zoom)  
FCI (Frequency Compound Imaging)  
B steer  
ExFOV (Extended Field of View)  
HR Flow™ (High Resolution Flow)  
Raw data processing  
iScanHelper  
1 active probe port  
Hard drive: 256 GB SSD or 1TB HDD  
4-USB  
HDMI  
iStorage  
MedTouch  
MedSight  
Net Storage  
Built-in Battery  
Power adapter  
Traveling case  
Multilingual controls overlay

### Optional Features

Free Xros M™  
Free Xros CM™  
Tissue Doppler Imaging  
TDI QA  
Continuous Wave Doppler

Low MI Contrast (Myocardium Contrast Imaging)  
Contrast Imaging QA (Quantitative Analysis)  
LVO (Left Ventricular Opacification)  
R-VQS (RF-Data based Quantitative Analysis on Vessel Stiffness)  
Strain Elastography  
Stress Echo  
Tissue Tracking QA  
RIMT  
AutoEF  
Auto DFR  
iWorks™ (Auto Workflow Protocol)  
iVocal  
McAfee  
DVR Module  
DICOM  
Clinical Measurement Package  
Mobile Trolley  
ECG module  
Internal Wi-Fi  
Dual-Probe extend module  
U-Bank (2 batteries)  
Barcode reader  
Footswitch  
iVocal Material package

### Language Support

Software: English  
Keyboard input: English  
Control panel overlay  
User manual

### Physical Specification

#### Dimensions and Weight

Width: 364±5 mm  
Depth: 322±5 mm  
Height: 44±3 mm  
Weight:  
About 3.0 kg (without battery)  
About 3.5 kg (with battery)

#### Monitor

15.6-inch high resolution color LED monitor  
Resolution: 1920 x1080  
Automatic brightness adjustment  
Screen Saver  
Open angle adjustable: 0 – 180°  
View angle (right/left): ≥170°

#### Handle

#### Probe Port

1 port connect to a transducer

#### Electrical Power

AC adapter Input:  
Voltage: 100 – 240 V AC  
Frequency: 50/60 Hz  
Power input: 2.0 – 1.0 A  
Battery: Lithium-Ion Battery Pack 14.4 V, 6600 mAh (single battery)



### Operating Environment

Ambient temperature: 0 – 40 °C  
Relative humidity: 20% – 85% (no condensation)  
Atmospheric pressure: 700 hPa – 1060 hPa

### Storage & Transportation Environment

Ambient temperature: -20 – 55°C  
Relative humidity: 20% – 95% (no condensation)  
Atmospheric pressure: 700hPa – 1060hPa

### User Interface

#### Control panel

Power/Battery Indicator  
Function Keys  
Ergonomic Soft Key Operation  
Backlit keys, ensuring accurate work in the dark room  
Programmable keys, available for userdefined functions  
Trackball, speed adjustment  
Key Brightness adjustment  
Integrated speakers, audio volume adjustment

#### Touchscreen

12.3-inch high sensitivity anti-glare color touchscreen  
Resolution: 1920x720  
Digital brightness and contrast adjustment through preset  
Viewing angle: ≥170 degrees  
Support touch screen gestures  
Support either hand writing or with gloves on

#### System Boot-up

Boot-up from complete shut-down: less than 26 sec  
Shut down: less than 30 sec

#### Comments

Supports text input and arrow  
Adjustable text size and arrow size and direction  
Supports home position  
Covers various application  
More than 800 comments items for versatile application  
User customizable

#### Bodymark

More than 232 bodymarks for versatile application

## Performance Specifications

### User Interface (continued)

#### Screen Information\* (presettable)

##### Common info:

Mindray logo  
Hospital name  
Exam date  
Exam time  
Acoustic power  
Mechanical index  
Tissue thermal index  
ID, Last name, First Name, Middle initial, Gender, Age  
Probe model  
ECG icon (when ECG connected)  
Operator  
TGC Curve  
Focus position  
Thumbnail  
Imaging parameters  
Help guidance  
Dynamic Trackball indices

\*Not all items are listed in this part, detail info please refer to user manual.

### Imaging Parameters

#### Overview

Digital beamformer  
Up to 1032192 channels  
64-beam forming

#### B-mode

Frame rate (max): 840 f/s  
A.Power: depend on probe  
TGC: 8 sliders  
Depth: 30 Levels  
Gain: 0 – 100, 1/step  
Steer: 5 Levels (available on linear transducers)  
FOV: On/off  
FOV Size: random adjustable  
FOV Position: random adjustable  
Image Quality: Pen/Gen/Res/HPen/HGen/HRes/HGen-FFR/HRes-FFR (depends on transducer)  
Persistence: 0 – 7, 1/step  
Dyn Ra.: 30 – 350 (depends on transducer)  
Gray Map: 1 – 8, 1/step  
Tint Map: Off, 1 – 8, 1/step  
ExFov: Off, 1 – 2 ((depends on transducer, extended FOV available on convex and linear transducers)  
iClear: Off, 1 – 7, 1/step  
iBeam: Off, 1 – 3, 1/step  
Line Density: L, M, H, UH  
L/R Flip: On/off  
U/D Flip: On/off  
Rotation: 0, 90°, 180°, 270°  
iTouch: On/off  
iTouch: -12 – 12, 3 db/step  
LGC: 8 point

Dual Live: On/off  
Auto Merge: On/off (available on linear transducers)  
H Scale: On/off  
Echo Boost: Off, 1 (available on phased transducers)  
Smooth: 0 – 6, 1/step  
TSI (Tissue Specific Imaging): General, Muscle, Fluid, Fat  
Zoom Value: 0.8 – 10  
HDScope: Off, 1 – 3, 1/step  
V1:1: On/off (available on linear transducers)  
iNeedle:  
B/iNeedle (on/off)  
Needle Dir: Auto, Left, Right

#### THI and PSH

Available on all types of transducer  
Patent PSH™ technology, obtains purer harmonic, better contrast resolution, higher SNR, exceptional high frequency harmonic  
iClear™ available  
Image quality: depends on transducers

#### M-mode

A.Power: depend on probe  
Gain: 0 – 100, 1/step  
Depth: same as B  
Speed: 25mm/s, 35mm/s, 50mm/s, 65mm/s, 100mm/s, 200mm/s  
Dynamic Range: 30 – 350 (depends on transducer)  
Gray Map: 1 – 8, 1/step  
Tint Map: Off, 1 – 8, 1/step  
Display format: V2:3, V3:2, H2:3, V3:1, FULL  
M Soften: 0 – 4, 1/step  
Edge Enhance: 0 – 3, 1/step  
Color M-mode available (convex and phased probe only)

#### Free Xros M

Speed: 25mm/s, 35mm/s, 50mm/s, 65mm/s, 100mm/s, 200mm/s  
Tint Map: Off, 1 – 8, 1/step  
Display Format: V2:3, V3:2, H2:3, V3:1  
Color Free Xros M available  
Gra Map: 1 – 8, 1/step  
Angle: adjustable  
Display: Cur./All; show A/B/C On/Off

#### Free Xros CM

Only available on TDI  
Speed: 25mm/s, 35mm/s, 50mm/s, 65mm/s, 100mm/s, 200mm/s  
Tint Map: Off, 1 – 8, 1/step  
Display Format: V2:3, V3:2, H2:3, V3:1  
Gray Map: 1 – 8, 1/step  
Angle: adjustable  
Edit, Undo, Delete function for curved line

#### Color Doppler Imaging

Frame rate (max): 331 f/s  
PRF: 0.1 kHz – 14.3 kHz

Velocity: 1.0 cm/s – 148.9 cm/s  
HR Flow™: High Resolution Flow provides better image quality and flow sensitivity  
A.power: same as B  
Gain: 0 – 100, 2/step  
Baseline: -8 – 8, 1/step  
Scale: 30 levels  
Quick Steer (available on linear transducers)  
Steer (available on linear transducers)  
ROI size/position: adjustable  
ROI Center Depth: adjustable  
Img Quality: Color/3 levels; HRFlow/1 level  
Persistence: 0 – 6, 1/step  
Smooth: 0 – 6, 1/step  
Color Map: V0 – V10; VV0 – VV9  
Flow State: L, M, H  
Priority: 0% – 100%, 1%/step  
WF: 8 Levels  
Line Density: L, M, H, UH  
Dual Live: On/off  
Invert: On/off  
Auto Invert: On/off (available on linear transducers)  
B/C Align: On/off  
Velocity tag: On/off  
Packet Size: 0 – 3, 1/step  
iTouch: On/off  
Smart Track: On/off

#### Power Doppler Imaging

PRF: 0.1 kHz – 14.3 kHz  
HR Flow™: High Resolution Flow provides better image quality and sensitivity  
A.power: same as B  
Gain: 0 – 100, 2/step  
Quick Steer (available on linear transducers)  
Steer (available on linear transducers)  
Scale: 30 steps  
ROI size/position: adjustable  
ROI Center Depth: adjustable  
Img Quality: Power/3 levels; HRFlow/1 level  
Persistence: 0 – 6, 1/step  
Smooth: 0 – 6, 1/step  
Dynamic Range: 10 – 70, 5/step  
Flow State: L, M, H  
Color Map: P0 – P3; dP0 – dP3  
Priority: 0% – 100%, 1/step  
WF: 8 levels  
Line Density: L, M, H, UH  
Dual Live: On/off  
Invert: On/off  
B/C Align: same as Color  
Packet Size: 0 – 3, 1/step  
iTouch: On/off  
Smart Track: n/off  
Auto Invert: On/off

## Performance Specifications

### Imaging Parameters (continued)

#### PW/CW-Mode

PW velocity:	11 cm/s – 770.0 cm/s
CW velocity:	5 cm/s – 3850.0 cm/s
PW PRF:	0.7 kHz – 20 kHz
CW PRF:	0.3 kHz – 100 kHz
A.Power:	same as B
Gain:	0 – 100, 2/step
Baseline:	9 levels
Quick Steer (available on linear transducers)	
Steer (available on linear transducers)	
Scale:	30 levels
Audio:	0% – 100%, 2%/step
Angle:	-89 – 89, 1/step
SVD:	random adjustable
Image Quality:	3 levels
Speed:	25mm/s, 35mm/s, 50mm/s, 65mm/s, 100mm/s, 200mm/s
SV:	0.5 – 30mm (PW only)
SV position:	random adjustable
Dynamic range:	24 – 72, 2/step
Gray map:	1 – 10, 1/step
Tint Map:	Off, 1 – 8, 1/step
Display format:	V2:3, V3:2, H2:3, V3:1, FULL
Invert:	On/off
Auto Invert:	On/off (available on linear transducers)
WF (depend on probe)	
Quick Angle:	-60°, 0°, 60°
Duplex/Triplex:	On/off
HPRF:	On/off
iTouch:	On/off
T/F Res:	0 – 6, 1/step
Auto Calculate:	On/off
Auto Calc Cycle:	1 – 5, 1/step
Trace Sensitivity:	0 – 5, 1/step
Auto Calc Parameter	
Trace Smooth:	Off, 1 – 4, 1/step
Trace Area:	Above, Below, All
Auto Calc Loop	

#### Tissue Velocity/Energy Imaging

Available on phased array transducer	
Max frame rate:	937.0 f/s
PRF:	0.4 kHz – 14.3 kHz
Velocity:	5 cm/s – 144.7 cm/s
A.Power:	same as B
Gain:	0 – 100, 2/step
Baseline:	-8 – 8, 1/step (TVI only)
Scale:	30 levels
Image Quality:	2 levels
Persistence:	0 – 6, 1/step
Smooth:	0 – 6, 1/step
Dyn Ra.:	10 – 70, 5/step (TEI only)
Tissue State:	L, M, H
Color Map:	
TVI:	TVV1 – TVV10
TEI:	P0 – P3, dP0 – dP3
Priority:	0 – 100, 1%/step
WF:	8 levels

Line Density:	L, M, H, UH
Dual live:	On/off
Invert:	On/off
B/C Align:	On/off
Velocity tag:	On/off (TVI only)
Packet size:	0 – 3, 1/step

#### Tissue Velocity Doppler

Available on phased array transducer	
Scale:	30 levels
Velocity:	7.01 cm/s – 616.0 cm/s
PRF:	0.7 kHz – 20 kHz
A.power:	same as B
Gain:	0 – 100, 2/step
Baseline:	9 levels
Audio:	0 – 100%, 2%/step
Angle:	-89 – 89, 1/step
SVD:	random adjustable
Img Quality:	2 levels
Speed:	25mm/s, 35mm/s, 50mm/s, 65mm/s, 100mm/s, 200mm/s
SV size:	same as PW
Dyn Ra.:	24 – 72, 2/step
Gray Map:	1 – 10, 1/step
Tint map:	Off, 1 – 8, 1/step
Display Format:	V2:3, V3:2, H2:3, V3:1, FULL
Invert:	On/off
WF:	10 levels
Quick Angle:	-60°, 0, 60°
Duplex/triplex:	same as PW
T/F Res:	0 – 6, 1/step
iTouch:	On/off

#### Tissue Velocity Motion

A.power:	same as B
Smooth:	0 – 6, 1/step
Velocity tag:	on/off
Persistence:	0 – 6, 1/step
Img Quality:	2 levels
Tissue State:	L, M, H
Speed:	25mm/s, 35mm/s, 50mm/s, 65mm/s, 100mm/s, 200mm/s
Display format:	V2:3, V3:2, H2:3, V3:1, FULL
Color Map:	TVV1 – TVV10
Packet Size:	0 – 3, 1/step
Priority:	0% – 100%, 1%/step
WF:	8 levels

#### Stress Echo

Available on cardiac sector transducers	
14 factory protocols	
User-defined protocols	
ECG triggered acquisition, display, selection, comparison, evaluation and archiving of multiple cardiac loops during various stages of a stress echo examination	
ASE16 (with score 4-7), ASE 17 (with score 4-7)	
Customized stages:	up to 7 views per stage, and up to 12 stages per study
View:	standard views (PSLA, PSAX, A4C, A2C), and customized views

#### Image acquisition:

R-wave trigger	
Acquire mode: Manual ROI or full screen	
Ability to acquire frames or clips in B-mode, LVO	

#### Image selection:

Attach the images with view annotation label (PSLA, PSAX, A4C, A2C, and customized views)

#### Review:

Automatically adjust to the number of images user defined

#### Wall Motion Scoring:

ASE 16 (with score 4-7), or ASE 17 (with score 4-7)	
Graphical display of scoring (Normal, Hyperkinetic, Severely Hyperkinetic, Akinetic, Dyskinetic)	

#### LV volume measurement

Measurement of LV Volume in all phases of cardiac cycle

#### Report

Reporting for both Wall Motion Scoring and LV volume measurement

Position: 0 – 100%

Curve Disp: 0.0 – 15.0

#### LVO

Available Probe:	SP5-1N5
Dedicated left ventricle contrast imaging tool	

#### iBeam™

Spatial compound imaging	
3 angles maximum	
Available on convex and linear transducers	

#### iTouch™

Auto image optimization	
B-mode:	gain, TGC
Color:	gain
Power:	gain
PW:	gain, scale, PRF, WF
Contrast imaging:	gain, TGC

#### Echo Boost™

Only for cardiac exams	
Improve the homogeneity of cardiac images through the whole field of view	
Better contrast resolution of myocardium tissue layers	
Better noise control in cardiac chambers and muscles	

#### B steer

Only for linear transducers

#### ExFov

Extended field of view	
Available for all convex and linear transducers	

#### Zoom

Zoom:	Spot zoom (write zoom) up to 10x, Pan zoom (read zoom) 0.8x – 10x
iZoom:	convertible 3 steps; normal image, zoom standard area, zoom only image area

#### QSave

Quick save image parameter setting after image adjustment done	
Support Save, Create, Restore	

## Performance Specifications

### Imaging Parameters (continued)

#### TDI QA

Dedicated quantification tool for TDI velocity, strain, strain rate analysis  
 Ellipse ROI, Standard ROI  
 Up to 8 of ROI  
 Delete all  
 Delete current  
 ROI tracking: tracking ROI along with cardiac movement  
 Smooth: 1 – 7, 1/step  
 X scale: 1 – 5, 1/step  
 Std.Height: 1.5 – 50 mm  
 Std.Width: 1.5 – 50 mm  
 Std.Angle: -89 – 90 degrees  
 Export: export current data as CSV format file

#### Tissue Tracking QA

Available probes:  
 SP5-1Ns  
 P8-2s  
 P10-4s  
 Tissue tracking quantitative analysis  
 Mandatory ECG connection before TT QA cine acquisition  
 Six views for analysis: ALAX, A4C, A2C, PSAXB, PSAXM, PSAXAP  
 Reload: reload cine again for new study  
 Edit: modify trace points  
 Start tracking  
 Accept & compute: start tracking myocardium movement when user accept trace result  
 Display effect: 0/1; at 0, tracking in velocity vector arrow; at 1, tracking in dots  
 Trace method: 3 point or manual for ALAX, A4C, A2C; manual for PSAXB, PSAXM, PSAXAP  
 Bull's Eye: trace result in bull's eye model  
 Torsion: Torsion rate curve display  
 LGC: available  
 Valve's open and close time index: MVC, MVC', AVC, AVO, MVO  
 Data export: export data in CSV file  
 Cycle: ECG triggered cardiac cycle recognition for analysis  
 Auto play: stop, X1/10, X1/5, X1/4, X1/3, X1/2, X1, X2, X3  
 Thickness: 1 – 30 mm, 1 mm/step; adjust trace thickness  
 Track point: 20 – 40, 1/step  
 Parameter: Volume, Speed, Displace., L Strain, L Strain R, T Strain, T Strain R, Area, R Strain, R Strain R, C Strain, C Strain R, Rotation, Rot. R  
 Smooth: 0 – 4, 1/step

#### AutoEF

Adjust Frame  
 Diastole FR  
 Systole FR

Volume curve: On/off  
 Adjustment for the border of endocardium

#### Smart Track

Continuously track the flow and detect the best color box position and angle in real time scanning.  
 The Linear probes in carotid, Upper Ext A, Upper Ext V, Lower Ext A, Lower Ext V, EM Vascular exam modes support the Smart Track function.

#### RIMT (RF-Data based IMT)

Available in single/dual B carotid exam mode  
 Side: left/right  
 Calculation of 6 RIMT values, RIMT average value, SD and ROI W  
 Report operation:  
 Data deleting  
 RIMT trend graphic viewing  
 Preview

#### Auto DFR

Automatic diastolic function assessment tool, automatically detect diastolic parameters of PW and TDI PW, to calculate a series diastolic function index E, A, E/A, e', E/e' automatically.  
 MV E/A  
 MV E/E' Septal  
 MV E/E' Lateral

### Cine Review and Post Processing

#### Cine review

Available in all modes  
 Frame by frame manual cineloop review or auto playback with variable speed  
 Independent cine review in 2D Dual and Quad mode one by one  
 Maximum cine memory is up to 25492 frames or 263.3 s (depend on the mode)  
 Retrospective storage (online setting available, 1 – 120 s, or 1 – 120 cycles, presettable) and prospective storage (1 – 480 s, or 1 – 390 cycles, pre-settable)  
 Frame compare: compare different frames for one cine in dual format  
 Cine compare: compare two or more than two cines in dual or quad format  
 Jump to first and jump to last: one keystroke review the first or last frame  
 Start point and end point: selectable

#### Raw data processing

B-mode:  
 TGC  
 Gain  
 Dynamic range  
 Gray map  
 Tint map  
 iClear  
 L/R Flip  
 U/D Flip  
 Rotation  
 LGC  
 Dual Live  
 Auto Merge  
 H Scale  
 Echo Boost  
 Smooth  
 Zoom Value

V1:1  
 B/iNeedle

#### M-mode:

Gain  
 Speed  
 Dynamic Range  
 Gray Map  
 Tint Map  
 Display format  
 Edge Enhance

#### Color:

Gain  
 Baseline  
 Smooth  
 Color map  
 Dual Live  
 Invert  
 Priority  
 Velocity tag

#### PW:

Gain  
 Baseline  
 Audio  
 Angle  
 Speed  
 Dynamic range  
 Gray map  
 Tint Map  
 Display format  
 Invert  
 WF  
 Quick Angle  
 T/F Res

### Measurement/Analysis and Report

#### Report

Specific report template by application  
 Editable value in report  
 Images selectable  
 Anatomy information  
 User-defined report template  
 Selecting report modules  
 Adding/removing measurement items from the report  
 Changing report layout  
 Load/save comment  
 Viewing history reports  
 Preview and printing reports  
 Able to Export as PDF file  
 Set the calculation method for the final value in batch

### Exam Storage and Management

#### Exam storage

SSD:  
 256 GB, more than 157 GB internal hard drive reserved for patient data storage  
 Capable of storage up to approximately 173242 single frames (FRM format)

## Performance Specifications

### Exam Storage and Management (continued)

#### HDD:

1 TB, more than 866 GB internal hard drive reserved for patient data storage

Capable of storage up to approximately 3290097 single frames (FRM format)

#### Storage area:

Pre-settable:	image area, standard area, full-screen
Image area:	1430x810
Standard area:	1920x920
Full-screen:	1920x1080

#### Exam management

iStation™ workstation dedicated for patient exam management

Patient exam query/retrieve

Support review of current and past exam

New exam, Active exam, Continue exam functions, End exam are available

Support measurements and calculations on archived exam and images

Export image as BMP/JPG/TIFF/DCM/FRM format (FRM: system format)

Export cine as DCM/AVI/CIN/MP4 format (CIN: system format)

Support backup/send to USB devices, DVD-RW media

#### iWorks™

Auto workflow protocol

Templates are user configurable

Functions: pause, stop, replace, repeat, skip, insert single step, return and continue, steps in thumbnail

iWorks setup mode: B/Dual/B+Color/ B+PW/ B+Color+PW/B+CW/ B+Color+CW/ B+M

iWorks setup annotation: support up to 2 annotations, location and font size are configurable

iWorks setup bodymark: select existing library, and probe indicator is pre-settable

iWorks setup measurement: select existing measurement library

Template import and export are available

### Connectivity

#### Ethernet Network Connection

Cable connection

Wireless connection: Internal Wi-Fi (including EAP Network)

#### DICOM 3.0

##### DICOM Basic

Verify (SCU, SCP)

Print

Store

Storage Commitment

Media Exchange

##### DICOM Worklist

##### DICOM Query/Retrieve

#### DICOM Modality Performed Procedure Step - MPPS

DICOM Cardiac structure report

DICOM Vascular structure report

#### iStorage

Direct network storage tool between ultrasound system and personal computer

#### MedSight

An interactive app that lets you transfer clinical images straight from Mindray Ultrasound system to a smart device, such as mobile phone or tablet PC

Needs to be installed on mobile terminal

Transfer images or clips from system to mobile terminal through WiFi

Support both iOS (7.0 and above) and Android (4.0 and above) system

For iOS powered smart device: DICOM is mandatory

For Android powered smart device: DICOM not necessary

#### MedTouch

Connect Ultrasound machine to smart devices based on Android and iOS system, such as tablet PC or mobile phone. Remote control of Ultrasound machine and tutorial software iScanHelper study on smart devices

Support Android and iOS powered smart devices

Android 4.0 and above

iOS 7.0 and above

DICOM not necessary

#### Net Storage

Support sending images or exams to the shared directory of your PC server.

#### iStorage

Data transfer

#### Security

Anti-Virus: McAfee and Windows Defender

VPN

### Transducers

#### Curved array

#### Linear array

#### L12-3RCs

Application: Abdomen, Pediatric, Small Parts, Musculo-skeletal, Vascular, Thoracic/Pleural

Bandwidth: 3.0 – 12.8 MHz

Depth: 1.5 – 35.0 cm

Number of Elements: 192

Field of View (max): 3.80 cm

Steered Angle: ±12°, ±6°, 0 (B steer); -30° – 30° (Color/PW steer)

Physical Footprint: 55.6x22 mm

Footprint: 43.5x8.2 mm

B-mode Frequencies: 3.0 – 8.3, 4.4 – 10.2, 5.6 – 12.8 MHz

Harmonic Frequencies: 8.0, 10.0, 12.0 MHz

Color Frequencies: 4.4, 5.0, 7.2, 7.2 (HR Flow) MHz

Doppler Frequencies: 4.2, 5.0, 7.2 MHz

Biopsy Guide: NGB-043 available, multiangle, reusable; NGB-044 available, multidepth, reusable

#### L9-3s

Application: Abdomen, Obstetrics, Pediatric, Small Parts, Musculo-skeletal, Vascular, Nerve

Bandwidth: 2.5-9.0 MHz

Depth: 1.5-35.0 cm

Number of Elements: 192

Field of View (max): 4.38 cm

Steered Angle: ±12°, ±6°, 0 (B steer); -10° – 10° (Color/PW steer)

Physical Footprint: 62x22 mm

Footprint: 48x11 mm

B-mode Frequencies: 2.5-7.0, 3.4-8.2, 3.6-9.0 MHz

Harmonic Frequencies: 5.0, 6.0, 7.0 MHz

Color Frequencies: 3.0, 3.6, 5.0, 7.3 (HR Flow) MHz

PW Frequencies: 3.0, 3.6, 5.0 MHz

Biopsy Guide: NGB-034, available, multiangle, reusable

#### Phased array

##### SP5-1Ns

Application: Abdomen, Gynecology, Obstetrics, Cardiac, Pediatric, Vascular, Thoracic/Pleural, Cephalic

Bandwidth: 1.5 – 4.5 MHz

Depth: 2.0 – 38.0 cm

Number of Elements: 64

Field of View (max): 90°

Physical Footprint: 38.2x30.4 mm

Footprint: 24x15.4 mm

B-mode Frequencies: 1.5 – 2.5, 2.5 – 3.5, 3.5 – 4.5 MHz

Harmonic Frequencies: 3.4, 3.8, 3.8, 4.2, 4.2 MHz

Color Frequencies: 2.0, 2.3, 2.5, 2.5 (HR Flow) MHz; TDI: 3.0, 3.8 MHz

PW Frequencies: 2.0, 2.3, 2.5 MHz; TDI: 2.5, 4.0 MHz

CW Frequency: 2.0 MHz

Biopsy Guide: NGB-011, available, multi-angle, reusable

##### P8-2s

Application: Abdomen, Pediatric, Cardiac, Cephalic

Bandwidth: 2.3 – 8.0 MHz

Depth: 2.0 – 38.0 cm

Number of Elements: 96

Field of View (max): 90°

Physical Footprint: 30.5x23.2 mm

Footprint: 19.5x11 mm

B-mode Frequencies: 2.3 – 5.4, 2.8 – 7.4, 4.2 – 8.0 MHz

Harmonic Frequencies: 5.0, 6.0, 7.0 MHz

Color Frequencies: 2.7, 3.3, 4.0, 2.5 (HR Flow) MHz; TDI: 3.0, 3.8 MHz

PW Frequencies: 2.7, 3.3, 4.0 MHz; TDI: 5.0, 6.0 MHz

CW Frequency: 2.5 MHz

Biopsy Guide: not available



## Performance Specifications

### Transducers (continued)

#### P10-4s

Application:	Abdomen, Pediatric, Cardiac, Nerve, Cephalic
Bandwidth:	3.0-11.4 MHz
Depth:	2.0-38.0 cm
Number of Elements:	128
Field of View (max):	90°
Physical Footprint:	15.1×10.2 mm
Footprint:	15×9.1 mm
B-mode Frequencies:	3.0-6.8, 3.8-10.2, 4.6-11.4 MHz
Harmonic Frequencies:	7.5, 8.0, 9.0 MHz
Color Frequencies:	4.0, 5.0, 5.7, 2.5 (HR Flow) MHz; TDI: 3.0, 3.8 MHz
PW Frequencies:	4.0, 5.0 5.7 MHz; TDI: 5.0, 5.7 MHz
CW Frequency:	5.0 MHz
Biopsy Guide:	not available

#### P7-3Ts

Application:	Cardiac
Bandwidth:	2.3-7.2 MHz
Depth:	2.0-38.0 cm
Number of Elements:	64
Field of View (max):	90°
Physical Footprint:	14×12 mm
Footprint:	12.2×12.2 mm
B-mode Frequencies:	2.3-5.4, 2.8-6.4, 3.3-7.2 MHz
Harmonic Frequencies:	6.0, 6.5, 7.0 MHz
Color Frequencies:	2.7, 3.3, 4.0, 2.5 (HR Flow) MHz; TDI: 3.0, 3.8 MHz
PW Frequencies:	2.7, 3.3, 4.0, 2.5 (HR Flow) MHz; TDI: 4.0 MHz
CW Frequency:	2.5 MHz
Biopsy Guide:	not available

#### Pencil

#### CW2s

Application:	Cardiac, Cephalic, Pediatric
Number of Elements:	2
Physical Footprint:	18.4×17.8 mm
Footprint:	18.4×17.8 mm
CW Frequency:	2.0 MHz
Biopsy Guide:	not available

### 10 Peripheral Devices and Accessories

#### Black/white digital video printer

MITSUBISHI P95DW-N

#### Black/white analog video printer

SONY UP-X898MD

#### Color digital video printer

SONY UP-D25MD

### Footswitch

USB port:	971-SWNOM (2-pedal/3-pedal)
USB port:	FS-81-SP (1-pedal)
Support User-definable functions (Freeze, Save, Print)	

### Built-in DVR

Built-in digital video recorder, save space and is a useful tool for education and memory  
Max storage length each time: 60 min

### Built-in Battery for Main Unit

Replaceable and rechargeable lithium battery  
Empty battery recharged to full in 4h  
Continuous work time: about 1.5 hour in B mode

### Mobile Trolley

#### MT3

Power supply module  
Dimensions (WxD): about 519 mm × 578 mm  
Platform height: 887 – 1207 mm; adjustable  
Weight:

- Without retractable cable and probe extend module: about 28.8 kg
- With retractable cable and without probe extend module: about 32.5 kg
- Without retractable cable and with probe extend module: about 30.9 kg
- With retractable cable and probe extend module: about 34.6 kg

#### Probe holders

- Auxiliary output cable
- Probe extend module
- Cover grounding cable
- Printer bracket

#### MT2

Dimensions (WxD): about 515 mm × 505 mm  
Platform height: 885 mm, 973 mm; 2 levels  
Weight:

- Without printer bracket and probe extend module: about 15 kg
- With printer bracket and probe extend module: about 18.8 kg
- Probe holders
- Probe extend module
- Printer bracket

### Barcode reader

- 1-D barcode reader: SYMBOL LS2208
- 2-D barcode reader: SYMBOL DS4308

### U-Bank

U-Bank with 2 batteries Weight: 1.95 kg

### Footswitch

USB port:	FS-81-SP-2(single pedal), 971-SWNOM (2/3-pedal)
Support User-definable functions (Freeze, Save, Print)	

### ECCG

6-pin, AHA/IEC, for 3-lead wires	
ECCG wave display:	On/off
ECCG source:	Lead/External
Position:	0 – 100%, 5%/step
Trig mode:	off/single/dual/timer
Gain:	0 – 30, 1/step
Sweep speed:	6 steps
Invert:	On/off

### Built-in Wireless adapter

Encryption:	WPA, WPA2
Max transfer speed:	300 Mbps
Protocols:	IEEE 802.11 ac/a/b/g/n
Frequency:	2.4G/5G

### System Inputs and Outputs

#### I/O Port

USB 3.0:	4 ports
Ethernet:	1 port
HDMI:	1 port
S-Video:	1 port

#### ECCG module

ECCG port:	1
------------	---

#### Probe Extend module

Probe port:	2
-------------	---

### Safety and Conformance

#### Quality standards

- ISO 9001
- ISO 13485

#### Design standards

- EN 60601-1 and IEC 60601-1
- EN 60601-1-2 and IEC 60601-1-2
- EN 60601-1-6 and IEC 60601-1-6
- EN 60601-2-37 and IEC60601-2-37
- EN 62304 and IEC 62304
- EN 62366 and IEC 62366
- EN ISO 17664 and ISO 17664

#### NOTICE:

Not all features or specifications described in this document may be available in all probes and/or modes. Mindray reserves the right to make changes in specifications and features shown herein, or discontinue the product at any time without notice or obligation.

Contact Mindray Representative for the most current information.