

P-Series Portable Ultrasound: Probes

The Sigma P5 Ultrasound offers superior imaging for multiple applications in Point Of Care, Musculoskeletal, Cardiovascular, and general imaging. The Lambda P9 Ultrasound is a fully functional cardiovascular echo system.


High Frequency Linear Probe (15L5V)

MSK, Vascular, and Small Parts applications
High definition, broadband probe.
Range: 15-5MHz
Center Frequency: 10 MHz
192 Elements



Phased array probe (2P4V)

Cardiac applications
Center frequency: 3 MHz
FOV (max): 90°




Linear array probe (12L5V)

MSK, Vascular and Small Parts applications
Center frequency: 7.5MHz
FOV: 38mm



Convex array probe (2C5V)

Abdominal, OB/GYN, and MSK applications
Center frequency: 3.5 MHz
FOV (max): 60°, Convex radius 60mm



Endo-Cavity (4EC9V)

OB/GYN and Prostate applications
Center frequency: 6.5MHz
FOV: 150.4°



Pediatric heart probe

Pediatric cardiac applications
Center frequency: 5.5MHz
FOV: 90°



Introducing the new wave of ultrasound



P-Series Ultrasounds