

E1 Digital B/W Portable Ultrasound System

Technical Specifications



SonoScape

SonoScape Medical Corp

1. General Specification

The system adopts the advanced ultrasonic Doppler technologies, including the Full Digital Super-wide Band Beam Former, Digital Dynamic Focusing, Variable Aperture and Dynamic Tracing, Wide Band Dynamic Range, and Multi-Beam Parallel Processing, etc. Users can perform the system with the minimum requirement of training or guidance. This system has been designed to comply with applicable international standards and regulations, ensuring the safety and availability of this product.

This system is based on the computer technology and Linux operation system, which make the system more flexible and stable. System maintenance and function update can be completed by software updating, through which would promote product value and keep the technological advancement.

2. Physical Specification

- Size: 378mm×352mm×114mm (W×H×D)
- Weight: Approx. 6.5kg (at most, including battery)
Approx. 6.1kg (at most, without battery)
- Monitor: 15.6'' widescreen and high-resolution color LCD monitor, LED backlight, anti-flickering and 45° angle adjustable
- Probe port: one (two ports can be equipped by order)

3. Advanced Technologies

- Digital front-end technology
- Multi-beam processing technology
- Spatial compound imaging

- μScan speckle reduction technology
- Tissue harmonic imaging
- Pure inversion harmonic imaging
- Graphic diagnosis icon

4. Standard Configurations

- μScan speckle reduction technology
- 5-band adjustable frequency in B mode
- Tissue specific imaging (TSI)
- THI mode
- PHI mode
- Quad beams
- 2D steer
- Spatial compound imaging
- Biopsy guide
- Spot zoom/ read zoom
- B/M/PW/CW auto optimization
- Basic measurement package
- Obstetrics measurement package
- Gynecology measurement package
- Cardiology measurement package
- Abdomen measurement package
- Vascular measurement package
- Urology measurement package
- Small parts measurement package
- Pediatrics measurement package
- TEI index
- Auto IMT measurement
- Standby mode
- Full screen zoom
- Show gallery
- DICOM3.0

5. Optional Configurations

- Needle visualization enhancement (linear)
- Trapezoidal imaging (convex)
- Panoramic imaging (2D)
- PW mode
- B+PW simult mode

- LGC
- PW auto trace

6. Optional Accessories

- ECG
- WiFi
- Biopsy bracket
- Color ink-jet printer
- B/W video printer
- Foot switch
- Trolley
- Backpack (with draw-bar)
- DVD record
- Remote control
- Barcode scanner
- Hard disk (1T)
- SSD hard disk (250G)
- Standard battery (90mins)
- Large capacity battery (120mins)

7. Scan Methods

- Electronic curved linear array scan
- Electronic linear array scan

8. Applications

- Abdominal
- Vascular
- Cardiac
- Gyn/OB
- Urology
- Musculo-skeletal
- Small parts
- Pediatric
- Cephalic
- Fetal
- Trans-rectal
- Trans-vaginal

9. Imaging Modes

- B mode
- M mode
- THI mode
- PHI mode
- PW mode
- B+PW simult mode
- Panoramic imaging (2D)

10. Display Formats (B Mode Includes Trapezoidal Imaging)

- B + B
- 4B
- B + M
- B + PW

11. System Parameters

- Frame rate: up to 60fps (probe dependent)
- Grayscale Level: 256

12. B Mode

- Gain: 0 - 255 adjustable, 5 per step
- Scan depth: 40cm
- Image zoom, showing magnification (0.8 - 10 times)
- TGC: 8 levels slider controls
- Image rotation: 0°, 90°, 180°, 270° selectable
- Steer: 0, ±2°, ±4°, ±6° selectable
- Harmonic imaging: Off, PHI, THI selectable
- Image inversion: Left and Right, Up and Down
- Compound imaging: Off, 1, 2, 3, 4 adjustable
- Focus: focus position and span adjustable, 17 levels adjustable. 1 represents single focus, 2 - 17 represents the length control of focus area.
- Frequency: 5 bands adjustable
- Chroma: Off and 12 types selectable
- μScan: Off, 1, 2, 3, 4, 5 selectable
- Line density: Low, Med, High adjustable (non-high

density)

- Persistence: Off, Low, Med, High, Max selectable
- Dynamic range: 20 - 200 selectable
- Grayscale curve: 16 selectable
- Sector width: adjustable
- Power: 1 - 100% adjustable, 5% per step
- RSI: adipose, muscle, fluid tissue and normal tissue
- Trapezoid imaging: On/Off
- Auto optimization function

13. M Mode

- Chroma: 13 types adjustable
- Display format: Full, H1/1, V1/2, V1/1, V2/1
- Scan speed: 5 levels adjustable
- Power: 1% - 100% adjustable

14. Spectral Doppler Mode

- Doppler methods: PW (pulsed wave) Doppler
- Mode: PW inactivated mode 1, inactivated mode 2, and activated mode.
- PW simult: On/Off (adjustable in inactivated mode 2 or activated mode)
- Sample volume and position for PW Doppler: 0.5 - 24.0mm adjustable
- Spectrum inversion: achievable
- Angle correction: 0°, 60°, -60° adjustable
- θ angle correction: -88° to 88°, 2° per step
- Spectral real-time trace: achievable
- Baseline shift: 9 steps selectable
- Frequency range: 3 steps selectable
- Wall filter: Min, Low, Med, High, Max adjustable
- PRF: 1 - 25KHz (PW) (probe dependent)
Max velocity range: ± 0.0576 to ± 6.14 m/s (PW) (probe dependent, sampling angle range: $\pm 60^\circ$)
- Scan speed: 5 levels adjustable
- Doppler Chroma: 13 types selectable
- Display format: FULL, H1/1, V1/2, V1/1, V2/1
- Steer angle: 0, $\pm 8^\circ$, $\pm 12^\circ$, $\pm 16^\circ$ adjustable

15. Integrated Data Management System

- Hard disk memory capacity: 500G (larger capacity optional)
- USB ports: three

16. Image Storage and Playback

- Cine playback: up to 10000 frames in B mode
- Static and dynamic image storage in B mode (including dual-split display and quad-split display)
- The stored images can be viewed directly on PC
- Clipboard function
- Doppler cine playback: Speed is adjustable; Sound can be played back

17. DICOM Network Communication

- Storage: directly transmits images with patient information to a DICOM file server
- Medical digital images and communication DICOM 3.0 interface

18. Preset Function

- Users can customize the presets based on different probes and examination parts to optimize imaging parameters
- Users can rearrange the presets.
- Users can import or export the presets. (not including default preset)

19. Patient Data Management

- Name, ID, Gender, Date of Birth, Height, Weight, LMP, EDD and GA can be input
- Patient reports and images are archived by patient examination
- Reports and images can be previewed
- Preview size can be set to 1×1 or 2×2

- Previewed file can be selected, deleted, printed or DICOM sent, DICOM printed or exported.
- Data can be exported to USB drive in BMP, JPG, TIF, AVI, WMV, PDF, TXT or HTML format.

20. Annotation and Body Mark Setting

- Body mark can be classified by specific examination. Body marks of the examination types such as abdomen, small part, urology, breast, gynecology, obstetrics, vascular, cardiology and MSK are provided.
- Annotation can be selected and input in the library. Annotations of the examination types such as abdomen, urology, small part, breast, gynecology, obstetrics, vascular, cardiology and MSK can be preset.

21. Safety Standard

- Comply with IEC 60601-1, Class I BF,
- Comply with IEC 60601-1-2, Group 1, Class B
- Comply with IEC 60601-2-37

22. Environmental Requirements

- Operation environment
 - Temperature: 0°C to +40°C
 - Relative humidity: 30% - 85% (no condensation)
 - Atmospheric pressure: 700 - 1060hPa
- Transportation and storage environment
 - Temperature: -20°C to +55°C
 - Relative humidity: 20% - 90% (no condensation)
 - Atmospheric pressure: 700 - 1060hPa
- Power supply
 - AC/DC Adaptor:
 - Input: 100-240V a.c., 1.5-0.75A, 50-60Hz;
 - Output: 19V d.c., 4.74A
 - Main unit: 19V d.c., 4.74A

23. Optional Probes

- Linear probes
 - L741
 - Frequency: 4.0 - 16.0MHz
 - Array element: 128
 - 10I2
 - Frequency: 4.0 - 16.0MHz
 - Array element: 96
 - L746
 - Frequency: 4.0 - 12.0MHz
 - Array element: 96
- Convex probes
 - 6V1
 - Frequency: 3.0 - 15.0MHz
 - Array element: 128
 - C361
 - Frequency: 2.0 - 6.0MHz
 - Array element: 96
 - C613
 - Frequency: 4.0 - 13.0MHz
 - Array element: 128
 - EC9-5
 - Frequency: 3.0 - 15.0MHz
 - Array element: 128

NOTE:

- The specifications of this system may change without any prior notification.
- Some products or features may not be available in some countries.
- Please contact your local SonoScape sales representative for more information.

Service Information:

Address: 4/F, 5/F, 8/F, 9/F & 10/F, Yizhe Building, Yuquan Road, Nanshan, Shenzhen, 518051, Guangdong, China

Tel: +86 755 2672 2890

Fax: +86 755 2672 2850

Email: Service@sonoscape.net

[Http://www.sonoscape.com](http://www.sonoscape.com)