



SiMPLi Series

SP-L01 Point-of-Care USB Linear Array Transducer

Ready for user-friendly ultrasound at the patient bedside? Simply connect the USB probe cable directly to your Windows laptop or tablet and you instantly have a lightweight, handheld and portable ultrasound imaging system for all of your point-of-care imaging needs.



Specifications

| | |
|---------------------|---|
| Depth Range | 0.1 cm - 10 cm |
| Focus | Dynamic over full depth range |
| Bandwidth | 5 - 10 MHz |
| Features | 2D, linear array, color flow imaging, lightweight transducer with flexible cable |
| Indications | Peripheral Vascular, Small Parts, Pediatrics, MSK, Neuromuscular, Pain Management, Acute Care |
| Weight | 80 grams |
| Cable | USB 2.0, 2 meter (5 meter optional) |
| System Requirements | OS: Windows 7, 8, 10 Processor/RAM: 2 GHz/4 GB |

Lightweight. Reliable. Flexible.

For more information regarding the SP-L01 or our other transducers, contact us at interson.com/contact-us or call us at **(925) 462-4948**.



Brachial Access



Thyroid/Carotid/Jugular



Carotid/Intima Media

