

M8 Elite Ultrasound System

Transducer Family



C11-3s (Micro-Convex)

Applications: Neurology, Pediatric Abdominal, Pediatric Cardiac, Vascular
FOV: 100°
Convex radius: 15 mm
Biopsy guide: NGB-018, multi-angle, reusable



L12-4s (Linear)

Applications: Musculoskeletal, Nerve, Pediatrics, Small Parts, Superficial, Vascular
FOV: 38 mm
Biopsy guide: NGB-007, multi-angle, reusable



C5-1s (Convex)

Applications: Adult Abdominal, Musculoskeletal, OB/GYN, Pediatrics, Vascular
FOV: 61°
Convex radius: 60 mm
Biopsy kit: NGB-022, multi-angle, reusable



L14-6Ns (Linear)

Applications: Musculoskeletal, Nerve, Neurology, Pediatrics, Small Parts, Superficial, Vascular
FOV: 38 mm
Biopsy guide: NGB-007, multi-angle, reusable



C6-2Gs (Convex)

Applications: OB, Abdominal, Pediatrics, Interventional
FOV: 92°
Convex radius: 20 mm
Biopsy guide: NGB-024, multi-angle, reusable



L16-4Hs (Linear)

Applications: Intra-operative, Musculoskeletal, Nerve, Superficial, Vascular Surgery
FOV: 25 mm
Biopsy guide: NA



D7-2s (Volume Convex)

Applications: Gynecology, Obstetrics, Abdomen
Bandwidth: 2.6 - 8.2 MHz
ExFOV: 110°
Radius: 40 mm
Needle guide: NA



L14-6Ws (Linear)

Applications: Breast, Musculoskeletal, Nerve, Small Parts
FOV: 50 mm
Biopsy guide: NGB-007, multi-angle, reusable



L9-3s (Linear)

Applications: Musculoskeletal, Vascular, Nerve
FOV: 43.8 mm
Biopsy guide: NGB-034, multi-angle, reusable



V11-3Ws (Endocavity)

Applications: OB/GYN, Urology
FOV: 173°
Convex radius: 11 mm
Biopsy guide: NGB-004, single-angle, reusable



L20-5s (Linear)

Applications: Musculoskeletal, Nerve, Small Parts, Superficial
FOV: 29 mm
Biopsy guide: NA



P7-3Ts (TEE)

Applications: Adult Transesophageal Echo
FOV: 90°
Biopsy guide: NA



SP5-1s (Single Crystal Phased)

Applications: Adult Abdominal, Adult Cardiac, Adult Transcranial, Pediatric Cardiac
FOV: 90°
Biopsy guide: NGB-011, multi-angle, reusable



P8-3Ts (TEE)

Applications: Pediatric Transesophageal Echo
FOV: 90°
Biopsy guide: NA



P10-4s (Phased)

Applications: Neonatal, Pediatric Abdominal, Pediatric Cardiac
FOV: 90°
Biopsy guide: NA



CW5s (Non-imaging Transducer)

Applications: Pediatric Cardiac, Vascular
CW frequency: 5.0 MHz
Biopsy guide: NA



P8-2s (Phased)

Applications: Pediatric Cardiac, Pediatric Abdominal
Bandwidth: 2.6 - 8.2 MHz
FOV: 90°
Needle guide: NA



CW2s (Non-imaging Transducer)

Applications: Adult Cardiac, Adult Cranial, Pediatric Cardiac
CW frequency: 2.0 MHz
Biopsy guide: NA

Mindray North America

800 MacArthur Blvd., Mahwah, NJ 07430

Tel: 800.288.2121 Support: 877.913.9663 Fax: 800.926.4275 www.mindray.com

Mindray® is a registered trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All brands and product names are trademarks of their respective owners. ©2024 Mindray DS USA, Inc. Subject to change. 04/24 P/N: 0002-08-40752 Rev 1.0

mindray
 healthcare within reach